


Kit Youcircuit Plus[®]

ANESTHETIC BREATHING CIRCUIT
BREATHING FILTERS

Pediatric/
Neonatal



Anesthetic Breathing Circuit

- Breathing bag
- ISO standard conical connections
- Adult and pediatric/neonatal versions
- Sterilized by Ethylene Oxide
- Surgical grade paper packaging
- Latex free

**SIZES UP TO
3 METERS**



Pediatric/Neonatal HMEF Filter - Filters and humidifies simultaneously

- Transparent body for better visualization
- Extensible corrugated tube
- ISO standard conical connections, 22/15 mm
- Connection for capnography

99.999%

OF EFFECTIVENESS
IN VIRUSES AND
BACTERIA FILTERING

Filters and protects

- Barrier filter, efficiency against viruses and bacteria
- Electrostatic and hydrophobic technology
- Transparent body for better visualization
- Connection for capnography

| Product/Description | Code |
|--|--------------|
| Pediatric/Neonatal Breathing Circuit 1,8 m, 1.5 m Extension, 1 l Bag and Connector 22M/22M + 1 HMEF Pediatric Breathing Filter with Catheter Mount | YCP-P18R15B1 |
| Pediatric/Neonatal Breathing Circuit 2 m, 1.5 m Extension, 1 l Bag and Connector 22M/22M + 1 HMEF Pediatric Breathing Filter with Catheter Mount | YCP-P2R15B1 |
| Pediatric/Neonatal Breathing Circuit 3 m, 1.5 m Extension, 1 l Bag and Connector 22M/22M + 1 HMEF Pediatric Breathing Filter with Catheter Mount | YCP-P3R15B1 |