


2nd
GENERATION

Youcircuit Plus®

ANESTHETIC BREATHING CIRCUIT



New Material Composition:

- Greater Resistance
- Better Handling
- More Safety

The **YOU CIRCUIT PLUS® ANESTHETIC BREATHING CIRCUIT** is a link between the patient and mechanical ventilation equipment and/or anesthesia machines.

It is used for the transportation of anesthetic gases, oxygen or other gases, in patients submitted to mechanical ventilation, when the maintenance of breathing is necessary for a long period of time.

The **YOU CIRCUIT PLUS® ANESTHETIC BREATHING CIRCUIT** can reach up to 3 m in length, and it can be retracted (corrugated tube) and return to original length, allowing the medical staff to adjust the best size of the circuit to carry out the procedure.


SIZES UP TO

3 METERS

Anesthetic Breathing Circuit

- Breathing bag
- ISO standard conical connections
- Adult and pediatric/neonatal versions
- Sterilized by Ethylene Oxide
- Surgical grade paper packaging
- Latex free

YOUCIRCUIT PLUS® ADULT



YOUCIRCUIT PLUS® PEDIATRIC/NEONATAL

Product/Description	Code	Anvisa
Adult Breathing Circuit 1,8 m	YCA+A18	81484769008
Adult Breathing Circuit 2 m	YCA+A2	81484769008
Adult Breathing Circuit 3 m	YCA+A3	81484769008
Pediatric/Neonatal Breathing Circuit 1,8 m	YCA+P18	81484769008
Pediatric/Neonatal Breathing Circuit 2 m	YCA+P2	81484769008
Pediatric/Neonatal Breathing Circuit 3 m	YCA+P3	81484769008
Adult Breathing Circuit 1,8 m, 1.5 m Extension, 3 l Bag and Connector 22M/22M	YCA+A18R15B3	81484769007
Adult Breathing Circuit 2 m, 1.5 m Extension, 3 l Bag and Connector 22M/22M	YCA+A2R15B3	81484769007
Adult Breathing Circuit 3 m, 1.5 m Extension, 3 l Bag and Connector 22M/22M	YCA+A3R15B3	81484769007
Pediatric/Neonatal Breathing Circuit 1,8 m, 1.5 m Extension, 1 l Bag and Connector 22M/22M	YCA+P18R15B1	81484769007
Pediatric/Neonatal Breathing Circuit 2 m, 1.5 m Extension, 1 l Bag and Connector 22M/22M	YCA+P2R15B1	81484769007
Pediatric/Neonatal Breathing Circuit 3 m, 1.5 m Extension, 1 l Bag and Connector 22M/22M	YCA+P3R15B1	81484769007