



You drain[®]

SILICONE CYLINDRICAL DRAIN


The **YOU DRAIN[®] SILICONE DRAIN** produced by ScavMedical is used for drainage of fluids and emptying of air accumulated in the body cavity in patients in the postoperative period.

It provides efficient fluid drainage, avoiding the accumulation of potentially harmful fluids from wounds, thus reducing the risk of infections and reinterventions.




Designed with a solid center to prevent collapse and obstruction, with "soft edges" to provide greater patient comfort.

ARTICULATING TROCAR FOR
BETTER ADJUSTMENT



Increased drainage safety

- Produced 100% in silicone
- Efficient and continuous drainage
- Star drain
- Radiopaque filament in all its extension (tube and drain)
- Anti-reflux valve
- Presentation with and without trocar



Product/Description	Code	Anvisa
Silicone Cylindrical Drain in 10 Fr with Bendable Trocar and Reservoir	YDS-10FR-KIT	81484760005
Silicone Cylindrical Drain in 15 Fr with Bendable Trocar and Reservoir	YDS-15FR-KIT	81484760005
Silicone Cylindrical Drain in 19 Fr with Bendable Trocar and Reservoir	YDS-19FR-KIT	81484760005
Silicone Cylindrical Drain in 10 Fr and Reservoir	YDS-10FR-KIT-DR	81484769005
Silicone Cylindrical Drain in 15 Fr and Reservoir	YDS-15FR-KIT-DR	81484769005
Silicone Cylindrical Drain in 19 Fr and Reservoir	YDS-19FR-KIT-DR	81484769005
Silicone Cylindrical Drain in 24 Fr and Reservoir	YDS-24FR-KIT-DR	81484769005