

Youlaparo®

FILTER FOR LAPAROSCOPY




The **YOULAPARO® FILTER FOR LAPAROSCOPY** is a passive system developed by ScavMedical, as an electrostatic and hydrophobic barrier filter, for use in minimally invasive procedures, in order to protect the surgical team of the suction of contaminated gases, expelled from the cavity of the patient to the external environment.

RATCHET CLAMP



Greater security in filtration and odor adsorption

- Efficiency in filtering viruses and bacteria: 99.999%
- Adsorbs up to 80% of odors
- Luer-Lock connector
- Ratchet clamp for Youlaparo fixing
- Double packaging, surgical grade paper and plastic film
- Latex free



**WITH
ACTIVATED
COAL**

The **YOULAPARO® FILTER FOR LAPAROSCOPY** has a filtration efficiency of 99.999% against viruses and bacteria and an activated carbon membrane which has the function of odor adsorption.

Product/Description	Code	Anvisa
Filter for Laparoscopy	YLP	81484760014