

# Youthoraxport®


DISPOSABLE THORACIC TROCAR



The **YOUTHORAXPORT® DISPOSABLE THORACIC TROCAR** was developed to be used in thoracic endoscopic procedures when insufflation is not necessary and in order to be a portal for the insertion of instruments and protection against foreign materials that may penetrate the thoracic cavity.

## Thoracic Trocar

- Indicated for use in the thoracic cavity
- Available in 4 diameters
- Disposable
- Blunt tip obturator
- Threaded cannula
- Silicone cover
- Latex free



Product/Description	Code	Anvisa
Disposable Thoracic Trocar 6 mm	YTP-06x7	81484769003
Disposable Thoracic Trocar 11 mm	YTP-11x7	81484769003
Disposable Thoracic Trocar 13 mm	YTP-13x7	81484769003
Disposable Thoracic Trocar 15 mm	YTP-15x7	81484769003