


Kit Youcircuit Plus[®]

ANESTHETIC BREATHING CIRCUIT
BREATHING FILTERS

Adult



Anesthetic Breathing Circuit

- Breathing bag
- ISO standard conical connections
- Adult and pediatric/neonatal versions
- Sterilized by Ethylene Oxide
- Surgical grade paper packaging
- Latex free

**SIZES UP TO
3 METERS**



HMEF Filter - Filters and humidifies simultaneously

- Transparent body for better visualization
- Extensible corrugated tube
- ISO standard conical connections, 22/15 mm
- Connection for capnography



99.999%

OF EFFECTIVENESS
IN VIRUSES AND
BACTERIA FILTERING

Electrostatic Bacterial/Viral Barrier Filter - Filters and protects

- Barrier filter, efficiency against viruses and bacteria
- Electrostatic and hydrophobic technology
- Transparent body for better visualization
- Connection for capnography

Product/Description	Code
Adult Breathing Circuit 1,8 m, 1.5 m Extension, 3 l Bag and Connector 22M/22M + 1 HMEF Adult Breathing Filter with Catheter Mount	YCP-A18R15B3
Adult Breathing Circuit 2 m, 1.5 m Extension, 3 l Bag and Connector 22M/22M + 1 HMEF Adult Breathing Filter with Catheter Mount	YCP-A2R15B3
Adult Breathing Circuit 3 m, 1.5 m Extension, 3 l Bag and Connector 22M/22M + 1 HMEF Adult Breathing Filter with Catheter Mount	YCP-A3R15B3
Adult Breathing Circuit 1,8 m, 1.5 m Extension, 3 l Bag and Connector 22M/22M + 1 HMEF Adult Breathing Filter with Catheter Mount + 1 Electrostatic Bacterial/Viral Barrier Filter	YCP-A18R15B3-E
Adult Breathing Circuit 2 m, 1.5 m Extension, 3 l Bag and Connector 22M/22M + 1 HMEF Adult Breathing Filter with Catheter Mount + 1 Electrostatic Bacterial/Viral Barrier Filter	YCP-A2R15B3-E
Adult Breathing Circuit 3 m, 1.5 m Extension, 3 l Bag and Connector 22M/22M + 1 HMEF Adult Breathing Filter with Catheter Mount + 1 Electrostatic Bacterial/Viral Barrier Filter	YCP-A3R15B3-E